

## Experiment 1: Newton's Ring

**AIM:** To measure the radius of curvature of a plano convex lens using Newton's rings apparatus

**APPARATUS:** (1) Newton's rings apparatus consisting of

- a. Plano Convex lens
  - b. Optically flat glass plate
  - c. Beam splitter
  - d. T-type traveling microscope with scale with L.C. = 0.001 cm
- (2) Monochromatic source of light of known wavelength (ex. Sodium)
- (3) Reading lamp and reading lens

### **OBSERVATIONS:**

**Table 1.1: Calculation of the least count of the scale on microscope**

Smallest Division on the main scale	0.05 cm
Number of Divisions on vernier scale	10
L.C. of traveling microscope	0.001 cm

**Table (1.2) Diameters of Newton's rings**

Seq. no. of Dark ring(n)	Upper position (P), cm	Lower position (Q), cm	Diameter $D_n = P - Q$ cm	Square of diameter $D_n^2, \text{cm}^2$
1	2.512	2.401	0.111	0.012321
2	2.534	2.381	0.153	0.023409
3	2.54	2.349	0.191	0.036481
4	2.553	2.553	0.218	0.047524
5	2.567	2.321	0.246	0.060516

### **CALCULATIONS:**

Slope of the graph of  $D_n^2$  Vs  $n = 0.01236$

Wavelength of sodium source used in the experiment= 5890 Å

Radius of curvature of plano convex lens

$$= \frac{\mu(Dm^2 - Dn^2)}{4(m-n)\lambda} = \frac{1 * slope}{4 * \lambda} = \frac{1 * 0.01236}{4 * 5890 * 10^{-8}} = 52.46 \text{ cm}$$

Standard radius of curvature $R_s$ cm	Radius of curvature using Newton's rings $R_e$ , cm	%deviation = $\left  \frac{R_s - R_e}{R_s} \right  * 100\%$
50	52.46	4.92 %

## GRAPH:

Reet Sharma  
Div-6 (F<sub>2</sub>)  
106030.


SCALE:

X-axis  $\rightarrow 2\text{cm} = 0.01\text{m}$

Y-axis  $\rightarrow 2\text{cm} = 1$ .

### Experiment - 01 Newton's Ring.

Graph of  $D_m^2$  Vs  $n$ .



SLOPE:

$$\begin{aligned} &= \frac{y_2 - y_1}{x_2 - x_1} = \frac{D_m^2 - D_m^2}{m - n} \\ &= \frac{0.0605 - 0.0234}{5 - 2} \\ &= \frac{0.0371}{3} \\ &= 0.01236 \end{aligned}$$

$\therefore$  Slope is 0.01236.

### **MY UNDERSTANDING OF THE EXPERIMENT:**

Newton's rings appear as a series of concentric, alternating bright and dark rings centred at the point of contact between two surfaces. We also observed that fringe width of the concentric circles decreases as we move in the outward direction.

In short, we learned the conditions to get Newton's Ring and measured the radius of curvature of a plano convex lens using Newton's rings apparatus

-----