Experiment 04 Laser based experiment I: Beam divergence

Performed On: 18/01/2022

 $\underline{\mathbf{AIM}} \boldsymbol{:}$ To measure the peak power and beam divergence of a given laser beam

<u>APPARATUS</u>: He-Ne laser, Optical bench, Laser Beam Analyzer with sensor and micrometerscrew arrangement.

Observation table 4.1				Observation table 4.2		
Powers at different positions at a distance $dI = 50 \text{ cm}$			Powers at different positions at a distance $d2 = 100 \text{ cm}$			
Sr. No.	Current in LBA, I (µA)	Position of micrometer x (mm)	Sr. No.	Current in LBA, I (µA)	Position of micrometer x (mm)	
i	0	-0.1	i	0	-0.1	
2	0	-0.075	2	0	-0.075	
3	0	-0.050	3	0	-0.05	
4	0	-0.025	4	0.03	-0.029	
5	0.06	-0.015	5	0.04	-0.025	
6	0.1	-0.011	6	0.06	-0.01	
7	0.13	-0.008	7	0.08	-0.015	
8	0.17	0	8	0.06	-0.025	
9	0.13	0.008	9	0.08	0.01	
10	0.10	0.011	10	0.06	0.015	
11	0.06	0.015	11	0.04	0.025	
12	0	0.025	12	0.03	0.029	
13	0	0.050	13	0	0.050	
14	0	0.075	14	0	0.075	
15	0	0.08	15	0	0.08	
16	0	0.1	16	0	0.1	

Calculations:

Divergence =
$$\frac{(D_2 - D_1) mm}{(d_2 - d_1) cm}$$

 $\frac{(D_2 - D_1) cm}{(d_2 - d_1)} \times 10^{-1}$
= $0.0005 \ rad$
= $0.00056 \ rad$ * $\frac{180 \ deg}{3.14 \ rad}$
= $0.028 \ deg$
= $1.68 \ min$

Table 4.3 Results

Sr.	Physlcal quantity	Value	Unit
No.			
1	Peak current the laser beam (at $d1=50$ cm)	0.17	mW
2	Divergence of laser beam	1.68	Min





Figure 4.3 Calculation of (a) diameter DI at position d1 (b) diameter D2 at position d2

Model Graph-I for Expt. 4, Laser Beam Divergence



Model Graph-II for Expt. 4, Laser Beam Divergence

